


Mineral Processing Lab



Gravity separation by Wilfley shaking table



Davis tube testing



Key Features & Capabilities

- **Sample Preparation:** Jaw Crusher (120mm coal to -30mm, -30mm to -4mm), Cone Crusher, Rod Mill, Ring Mill, eight segment Rotary Sample Dividers of various sizes, AS 1210 compliant Pressure Filters
- **Gravity Separation:** Wilfley Shaking Table (size #13 standard)
- **Magnetic Separation:** Davis Tube Tester
- **Flotation:** Denver Flotation Cells of various sizes
- **Size Classification:** Sieve Shaker (both dry and wet sieving), Cyclosizer (sub-sieve range)



Cyclosizer (originally known as the Warman Cyclosizer)



Past Sponsors --- Lithium Australia, Venus Metals Corporation, MRIWA, Arcadium Lithium, King River Resources, Podium Minerals, LKAB (Sweden), FBI CRC



Contacts

- **Lead:** Dr. Rorie Gilligan
- **Contact:** Prof. Aleksandar Nikoloski
a.nikoloski@murdoch.edu.au